

## DTPM 16-PANEL POINT-OF-CARE DRUG CUP RELIABLE URINE DRUG SCREENING IN JUST 5 MINUTES

In today's fast-paced world, reliable and efficient drug screening tools are essential for testing facilities. The DTPM 16-Panel POC Cup stands out by providing accurate, user-friendly, and cost-effective urine drug screening results in just five minutes.

The DTPM 16-Panel POC Cup was thoughtfully designed with features that make it one of the best options for quick, accurate, and effective drug screening.

### INTEGRATED PROCEDURAL CONTROLS:

- ✓ Minimum and maximum fill lines help ensure the proper sample volume was collected.
- ✓ Equipped with a temperature strip, the DTPM 16-Panel POC Cup ensures that samples measure within the acceptable temperature range, adding an extra layer of accuracy to the screening process.
- ✓ Internal quality controls on each test strip confirm sufficient sample volume, adequate wicking, and correct procedural technique were used providing additional confidence in the test's accuracy.
- ✓ Sample adulteration is evaluated by four integrated Sample Validity Tests: Creatinine, Oxidants, pH, and Specific Gravity.



SCAN TO ORDER  
SKU: DTPM-16PV2



### THE DTPM 16-PANEL POC CUP INCLUDES THE FOLLOWING TESTS:

DRUGS OF ABUSE	CUT-OFF LEVEL
AMP = Amphetamine	1000 ng/mL
BAR = Barbiturates	300 ng/mL
BUP = Buprenorphine	10 ng/mL
BZO = Benzodiazepines	300 ng/mL
COC = Cocaine	300 ng/mL
ETG = Ethyl Glucuronide	500 ng/mL
FTY = Fentanyl	200 ng/mL
K2 = Synthetic Cannabinoids	50 ng/mL
MDMA = Ecstasy	500 ng/mL
MET = Methamphetamine	1000 ng/mL
MOP = Morphine	300 ng/mL
MTD = Methadone	300 ng/mL
OXY = Oxycodone	100 ng/mL
TCA = Tricyclic Antidepressants	1000 ng/mL
THC = Tetrahydrocannabinol	50 ng/mL
TRA = Tramadol	200 ng/mL
ADULTERANTS:	CR; OX; pH; SG

### QUICKLY TRAIN STAFF USING THE ONLINE TRAINING VIDEO WITH OPTIONAL CERTIFICATION QUIZ.



VISIT [DTPM.COM/POC-CUP](http://DTPM.COM/POC-CUP)

**DISCOUNTS AVAILABLE FOR LARGE QUANTITY ORDERS**

Call us today for a customized pricing quote!