

ELIMINATE URINE SPILLS, IMPROVE EFFICIENCY, REDUCE COSTS

WITH

**STERILE VACUUM
URINE CONTAINER**
WITH INTEGRATED FILLING PORT

TA990001
PACK OF 250



**VACUUM URINE
COLLECTION TUBE**

TA990002
PACK OF 100

TA990002CS
CASE OF 1800

EASY TO USE:

- 1 Collect urine in container.
- 2 Simply press the urine tube into the filling port – the tube will automatically fill.
- 3 Use tube for testing, storing, and transporting.

BENEFITS:

- ✓ Eliminate urine pour over processes.
- ✓ Reduce storage space requirements.
- ✓ Cleaner & more cost-effective than urine cups.
- ✓ Improved workflow reducing hands-on processes.



DTPM

SHOP ALL SUPPLIES AT
STORE.DTPM.COM

CONTACT SALES
256-845-1533