

EXPANDED CONFIRMATION TEST MENU FROM DTPM'S INDUSTRIAL LABORATORY

DTPM, Inc is a nationwide provider of laboratory products, services, and solutions. DTPM has partnered with laboratories since 1993 to provide a comprehensive array of instrumentation, custom assays, laboratory consumables, laboratory supplies, and reliable service.

DTPM announces the availability of new confirmatory testing for the following drugs using LC-MS/MS. Many of these drugs are prescribed by physicians but are often seen in substance use situations. The drugs and their respective cutoffs are listed below.

Drug Name	Cutoff Level
Fentanyl/Norfentanyl	2.0 ng/mL
Xylazine	50 ng/mL
Mitragynine (Kratom)	10 ng/mL
Carisoprodol (Soma® , Vanadom®)	50 ng/mL
Pregabalin (Lyrica®):	500 ng/mL
Gabapentin (Neurontin®)	500 ng/mL
Tianeptine & MC5 Metabolite	10 ng/mL
Cyclobenzaprine (Flexeril® , Amrix®)	50 ng/mL
Ritalinic Acid* (Ritalin® , Methylin®, Concerta®)	50 ng/mL

*This is the major metabolite to accompany our current Methylphenidate test

TEST REQUIREMENTS

Sample Type: Urine

Collection Device: Standard urine collection cup

Sample Storage: Refrigerate after collection prior to shipping to DTPM. Storage in extreme heat or uncontrolled temperatures for longer than seven days may result in decreased accuracy.

Sample Shipping: Ship sample and chain of custody documentation overnight to DTPM for best results. Biological shipping status must be used.

Rejection Criteria: Samples that have illegible, contradictory, or missing information on the chain of custody form may be rejected. Sample leakage may result in rejection.

Contact us today for more information at 256-845-1261.

Visit the DTPM Store at Store.DTPM.com for drugs of abuse screening reagents, urine collection cups, and general supplies such as gloves, lint-free wipes, transfer pipettes, etc.

Call us today for a customized pricing quote!