

## ELIMINATE URINE SPILLS, IMPROVE EFFICIENCY, AND REDUCE COSTS

*DTPM, Inc is a nationwide provider of laboratory products, services, and solutions. DTPM has partnered with laboratories since 1993 to provide a comprehensive array of instrumentation, custom assays, laboratory consumables, laboratory supplies, and reliable service.*

Reduce storage space requirements and eliminate sample pour-over processes with our Sterile Vacuum Urine Container with Integrated Filling Port + Vacuum Urine Collection Tube. These two products will improve your testing workflow as they are cleaner and more cost-effective than standard urine cups.

### Easy to Use:

- Collect urine in container.
- Simply press the urine tube into the filling port – the tube will automatically fill.
- Use tube for testing, storing, and transporting.

### Benefits:

- Eliminate urine pour over processes.
- Reduce storage space requirements.
- Cleaner & more cost-effective than urine cups.
- Improved workflow & happier lab personnel!



### Associated Part Numbers:

**Sterile Vacuum Urine Container with Integrated Filling Port**

TA990001: *Pack of 250*

**Vacuum Urine Collection Tube**

TA990002: *Pack of 100*

TA990002CS: *Case of 1800*

**Contact us today for more information at 256-845-1261.**

Visit the DTPM Store at [Store.DTPM.com](http://Store.DTPM.com) for drugs of abuse screening reagents, urine collection cups, and general supplies such as gloves, lint-free wipes, transfer pipettes, etc.

**Call us today for a customized pricing quote!**