

CLIA-WAIVED 13-PANEL POINT OF CARE CUPS

DTPM
TOXICOLOGY

THE DTPM SERVICE DIFFERENCE:

DTPM recently designed a 13-Panel CLIA-waived POC cup that is proven to perform superior and is more accurate than the leading alternate brand POC cup. The DOA 13-panel Test Cup incorporates a one-step approach into a total system for on-site urine drug testing. An accurate alternative to laboratory testing with results in 5 minutes, not 24-48 hours as required by laboratories, and is one of the most accurate tests available.

POC CUP FEATURES INCLUDE:

- No tipping or tilting – Our drug test cups are the easiest in the market!
- No dipping reagent strips or test card – A fully integrated test that provides solutions in minutes. No training, special skills, or help from a third party is required to administer the drug tests.
- Accuracy – With comparable accuracy to other brands and FDA approval, you are likely to have the correct results if the test is performed properly.

THE DTPM CLIA-WAIVED 13-PANEL POC CUP TEST FOR THE FOLLOWING:

DRUGS OF ABUSE:	CUT-OFF LEVEL
AMP = Amphetamine	1000ng/ml
BAR = Barbiturates	300ng/ml
BUP = Buprenorphine	10ng/ml
BZO = Benzodiazepines	300ng/ml
COC = Cocaine	300ng/ml
MDMA = Ecstasy	500ng/ml
MET = Methamphetamine	1000ng/ml
MOP = Morphine	300ng/ml
MTD = Methadone	300ng/ml
OXY = Oxycodone	100ng/ml
PCP = Phencyclidine	25ng/ml
TCA = Tricyclic	1000ng/ml
THC = Tetrahydrocannabinol	50ng/ml
ADULTERANTS:	CR; OX; pH; SG



SCAN TO VISIT
OUR DTPM STORE



Call us today for a customized pricing quote!