

QUICK AND CONVENIENT 16-PANEL POINT OF CARE CUPS

DTPM
TOXICOLOGY

THE DTPM SERVICE DIFFERENCE:

DTPM recently designed a new DTPM 16-Panel POC cup that is proven to perform superior and is more accurate than the leading alternate brand POC cup. The DOA 16-panel Test Cup incorporates a one-step approach into a total system for on-site urine drug testing. An accurate alternative to laboratory testing with results in 5 minutes, not 24-48 hours as required by laboratories, and is one of the most accurate tests available.

POC CUP FEATURES INCLUDE:

- No tipping or tilting – Our drug test cups are the easiest in the market!
- No dipping reagent strips or test card – A fully integrated test that provides solutions in minutes. No training, special skills, or help from a third party is required to administer the drug tests.
- Accuracy – With comparable accuracy to other brands and FDA approval, you are likely to have the correct results if the test is performed properly.

THE DTPM 16-PANEL POC CUP TEST FOR THE FOLLOWING:

DRUGS OF ABUSE: CUT-OFF LEVEL

| | |
|-----------------------------|----------------|
| AMP = Amphetamine | 1000ng/ml |
| BAR = Barbiturates | 300ng/ml |
| BUP = Buprenorphine | 10ng/ml |
| BZO = Benzodiazepines | 300ng/ml |
| COC = Cocaine | 300ng/ml |
| ETG = Ethyl Glucuronide | 500ng/ml |
| FTY = Fentanyl | 200ng/ml |
| K2 = Synthetic Cannabinoids | 50ng/ml |
| MDMA = Ecstasy | 500ng/ml |
| MET = Methamphetamine | 1000ng/ml |
| MOP = Morphine | 300ng/ml |
| MTD = Methadone | 300ng/ml |
| OXY = Oxycodone | 100ng/ml |
| TCA = Tricyclic | 1000ng/ml |
| THC = Tetrahydrocannabinol | 50ng/ml |
| TRA = Tramadol | 200ng/ml |
| ADULTERANTS: | CR; OX; pH; SG |



SCAN TO VISIT
OUR DTPM STORE



Call us today for a customized pricing quote!